Aquamantys® 6.0 Bipolar Sealer Aquamantys® 2.3 Bipolar Sealer

REF 23-112-1

REF 23-113-1

70-10-1414 Rev D
Eff. Date: May 2012

Follow Instructions for Use

Do Not Reuse

Do Not Resterilize

Humidity Limitation
Operate at 15% – 85%
Non-condensing Humidity

Sterilized by Ethylene Oxide

Does not contain
natural rubber latex

Contains Phthalates

Use by Date

Temperature Limitation
Operate at 50 °F to 86 °F
(10 °C to 30 °C)

Non-Pyrogenic

Medtronic

Manufacturer:

Medtronic Advanced Energy LLC

180 International Drive

Portsmouth, NH 03801 USA

Tel 866.777.9400

Fax 866.222.0900

www.medtronic.com/advancedenergy

Medtronic B.V.

Earl Bakkenstraat 10

6422 PJ Heerlen

The Netherlands

Medtronic Australasia Pty Ltd.

97 Waterloo Road

North Ryde NSW 2113 Australia


© Copyright 2011 Medtronic Advanced Energy LLC. All rights reserved. Printed in USA. Aquamantys®, the Aquamantys® logo, and Transcollation® are registered trademarks of Medtronic, LLC. Additional trademarks are the property of their respective owners.

ENGLISH


Device Description

The Aquamantys® Bipolar Sealer is a sterile, bipolar device. The device employs radio-frequency (RF) energy and saline irrigation for hemostatic sealing and coagulation. The device is equipped with a dual electrode tip. Saline and electrical lines exit the opposite end of the handpiece from the electrodes. The handpiece is equipped with an on-off button that simultaneously activates both RF and saline flow. A saline fluid delivery line is provided with the device, and includes a section of pump tubing and a drip chamber or spike. The three-pin electrical connector is designed to be plugged into the Aquamantys Pump Generator.

Figure 1. Side view of Aquamantys Bipolar Sealer (Length of cable and saline lines not to scale)

System Diagram: The device is to be used only with the Aquamantys Pump Generator.

Indications for Use

The Aquamantys Bipolar Sealer is a sterile, bipolar electrosurgical device intended to be used in conjunction with the Aquamantys Pump Generator for delivery of RF energy and saline for hemostatic sealing and coagulation of soft tissue and bone at the operative site. It is intended for, but not limited to orthopaedic, spine, thoracic, and open abdominal surgery. The device is not intended for contraceptive tubal coagulation (permanent female sterilization).

Warnings

- Special care should be taken when using the device in the proximity of neural tissue to avoid damage to nerves and similar sensitive structures.
- DO NOT TREAT the following areas at any time:
 - Dura
 - Nerve roots
 - Skin and skin edges
 - Intact nerves
 - Intact tendons and ligaments
 - Epidural veins
 - Bone surfaces that are intended to be fused
 - Vertebral end plates after discectomy
 - Subcutaneous tissue
 - Bone to be covered by implant
- The device is not intended for contraceptive tubal coagulation (permanent female sterilization).
- **Ensure that both electrodes are in contact with the tissue to be treated. Activation and saline flow occur simultaneously. Ensure that saline is flowing at the time of activation.** Activating device with tips pointing straight upward may result in inadequate saline flow to surface intended for treatment.
- **If saline flow stops during the electrosurgical procedure, stop using the device and attempt to resume saline flow. Ensure that the pump tubing segment has been loaded properly into the pump head located on the generator and that the saline bag is not empty. If unable to resume saline flow, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**
- This device is not intended to be bent, or used as a pry, or any other use not cleared by the FDA. Bending or using the device as a pry could cause part breakage.
- The device is provided as a sterile, non-pyrogenic, single-use device. **Do not resterilize or reuse this device.** Reprocessing (resterilizing or reusing) this device can result in occluded saline apertures, reducing or preventing the flow of saline.
- Be aware that all exposed metal on the electrodes is capable of treating tissue. Use caution to avoid inadvertent treatment of tissue and adjacent structures.
- Be aware that the device employs RF coupled with saline. This coupling effect may result in a deeper tissue effect than conventional RF and has the potential for hot saline run-off onto delicate structures.
- Protect delicate structures from hot saline run-off by utilization of suction or other protective measures.
- Use suction to avoid activating the device in a pool of saline. Activating in pooled saline may reduce the hemostatic effectiveness of the device.
- Use suction to minimize the potential for activation of the device in a pool of blood. Activating in a pool of blood may limit the hemostatic effectiveness of the device or increase the risk of an electrode becoming dogged by coagulated blood.
- The tip of the suction wand should not touch the electrodes, since this might interfere with the proper function of the device. However, the tip of the suction wand can be as close as 1 – 2 mm from the electrodes when the device is activated.
- Overuse or excessive application of this device may result in contraction, inflammation or necrosis of tissue.
- Use of this device on skin may result in incisional complications such as necrosis or desiccation of the skin.
- Inspect the device and cord for breaks, cracks, nicks, or other damage before use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team.
- Interference produced by the operation of the RF surgical equipment may adversely influence the operation of other electronic equipment.
- Failure of the RF surgical equipment could result in an unintended increase of output power.
- Use the lowest setting possible to achieve the desired tissue effect to avoid overtreatment which could result in swelling, fluid, seroma or unintended tissue necrosis.

Precautions

- Surgery should be performed by persons with adequate training and preparation. Personnel should fully understand the nature and use of RF before performing electrosurgical procedures to avoid the risks of shock and burn hazards to both the patient and the operator and damage to the instrumentation.
- The cable on the device should be positioned in a way to avoid contact with the patient or other cables.
- Consult the operating and user manuals for light sources, other electrosurgical units, and other ancillary devices for operating instructions, warnings, and cautions prior to their use in the same surgical field as Medtronic's device.
- It is recommended that physicians utilize pre-clinical training, review of pertinent literature and other appropriate educational tools before attempting newer surgical procedures, such as endoscopic, laparoscopic, or thoracoscopic procedures.
- Examine the shipping carton, packaging, sterile barrier, and device for any signs of transit damage. If there are any shortages, breakage, or apparent damage, do not use the device. Return the device to Medtronic Advanced Energy LLC and use a new device.
- This device does contain Phthalates.
- **Use this device only with the Aquamantys Pump Generator.** Read the warnings, precautions, and instructions provided with Aquamantys Pump Generator before using. Specific instructions are not included in this manual.
- Use the device with caution in the presence of pacemakers, as electrosurgical devices may cause interference with pacemakers or other active implants.
- Place any monitoring electrodes being used as far away as possible from the device to avoid electrical interference with monitoring equipment.
 - Avoid needle-monitoring electrodes.
 - Use monitoring systems incorporating high frequency current limiting devices.
- The depth of effect is deeper and increases with time if the electrodes are held stationary with less depth of effect if the electrodes are moved over tissue.
- High power settings result in deeper tissue effect than lower power settings.
- DO NOT use electrosurgery in the presence of flammable anesthetics or other flammable gases, near flammable fluids or objects, or in the presence of oxidizing agents, as fire could result.
- The patient should not come into contact with metal parts which are earthed or have an appreciable capacitance to earth (e.g. operating table supports, etc.)
- Skin to Skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.
- For surgical procedures where the RF current could flow through parts of the body having a relatively small cross sectional area, the use of bipolar techniques may be desirable in order to avoid unwanted tissue damage.
- Before using the device, confirm the following:
 - The cable from the device is connected to the Aquamantys Pump Generator.
 - All electrical connections are tight, clean, and dry.
 - All fluid connections are secure.
 - The generator is set at the desired power level.
 - The saline delivery tubing and device have been fully primed with sterile saline (0.9% NaCl) solution.

Setting Up the Aquamantys Disposable Bipolar Device

Warnings:

- **Electric Shock Hazard.** Ensure that the device plug is correctly connected and that no metal pins are exposed.

Precautions:

- Read the instructions, warnings, and precautions provided with the electrosurgical device and the Aquamantys Pump Generator before using.

STEP 1: Connect the Aquamantys Disposable Bipolar Device to the Aquamantys Pump Generator

1. Using aseptic technique, open the device package and deliver contents to the sterile field.
2. Using aseptic technique, pass the capped drip chamber or spike of the saline line off the sterile field, making sure to maintain an adequate length of saline line and electrical cable on the sterile field.
3. Connect the Aquamantys disposable bipolar device to the Aquamantys Pump Generator by inserting the plug of the device into the plug receptacle on the front panel of the pump generator.

STEP 2: Load the pump segment portion of the Aquamantys Disposable Bipolar Device into the pump head of the Aquamantys Pump Generator

Warning:

- Always close the pump head prior to priming or device activation. Always allow the pump head rotor to come to a complete stop prior to opening the pump head. Do not attempt to load or adjust the positioning of the pump segment of the Aquamantys disposable bipolar device in the pump head while the pump head rotor is turning. Fingers or loose clothing could be caught in the pump rollers.

The saline delivery tubing of the Aquamantys disposable bipolar device includes a special pump segment portion designed to operate with the pump head of the Aquamantys Pump Generator. The pump segment portion of the saline delivery tubing is located between a black tubing connector and a white tubing connector.

The pump head is located on the right side of the Aquamantys Pump Generator when looking at the unit from the front. It is best to position yourself facing the right side of the unit to load the pump segment portion of the Aquamantys disposable bipolar device into the pump head.

1. Use the black-tipped lever located on the right side of the pump head to open the pump head. Rotate the black-tipped lever 180° (degrees) counterclockwise from the right side of the pump head to the left side of the pump head. This action will raise the upper part of the pump head.
2. After locating the pump segment portion of the saline delivery tubing on the Aquamantys disposable bipolar device, place the pump segment portion of the saline delivery tubing into the pump head with the black tubing connector positioned to the left side of the pump head (i.e. closest to the front panel of the Aquamantys Pump Generator). The white tubing connector should then be positioned to the right side of the pump head.

Warning:

- Do not peel saline delivery segment apart from the electrical cable before placing the pump segment in the pump head. Peeling the tubing first increases the potential for loading the pump segment in the reversed position.
3. Use the black-tipped lever to close the pump head. Rotate the black-tipped lever 180° (degrees) clockwise from the left side of the pump head to the right side of the pump head. This action will lower the upper part of the pump head.
 4. At the locations where the tubing enters and exits the pump head, the upper (moving) part of the pump head includes black slotted tubing guides. Ensure that the pump segment portion of the saline delivery tubing is properly aligned in the pump head by inspecting where the tubing enters and exits the pump head. The pump segment must be centered in the guide slot of both tubing guides, with no pinching of the tubing.

STEP 3: Spike the saline bag

1. Hang a 500 ml or larger bag of sterile saline (0.9% NaCl) solution on the Aquamantys Cart I.V. pole or another I.V. support, which is in close proximity to the Aquamantys Pump Generator.
2. Remove the protective cover over the spike at the end of the Aquamantys disposable bipolar device saline delivery tubing.
3. Using aseptic technique, spike the bag of sterile saline (0.9% NaCl) solution.
4. If a drip chamber is present, squeeze the drip chamber once or twice to fill the drip chamber to a level of at least one third full.
5. Open the vent cap if the source of sterile saline (0.9% NaCl) solution is a non-vented glass bottle.

STEP 4: Prime the Aquamantys Disposable Bipolar Device

Press the START PRIME button on the front panel of the Aquamantys Pump Generator. This initiates priming of the Aquamantys device with saline. Refer to the Aquamantys Pump Generator User Guide.

Warning:

- Temporarily unused active electrodes should be stored in a location that is isolated from the patient. Always place the device into a holster when not in use to avoid unintended activation that can cause injury to the patient or surgical team.
- Always place the device into a holster or over a container to collect the saline that exits the electrodes as a result of the priming process. If excess saline is not collected, saline could drip on the patient, patient drapes, surgical instruments, or operating room surfaces.
- Lack of saline flow from both of the electrodes can result in a lack of tissue effect and may damage the electrodes during device activation. Use caution to avoid the following conditions that can result in lack of adequate saline flow from the device.
 1. Pump segment portion of the saline delivery tubing loaded improperly into the pump head:
 - a. In the wrong direction. The black tubing connector should be to the left side of the pump head (i.e. closest to the front panel of the Aquamantys Pump Generator.)
 - b. Pinched pump segment portion. Tubing not aligned in the guide slot of each tubing guide.
 - c. Upper part of pump head not completely lowered onto the pump segment portion of the saline delivery tubing. The upper part of the pump head must be completely lowered all the way down (black lever rotated all the way to the right), so that the pump head can properly interact with the pump segment portion of the saline delivery tubing.
 - d. Pump segment not loaded into the pump head at all.
 2. Priming not completed:
 - a. START PRIME button not pressed.
 - b. START PRIME button pushed before the saline bag was spiked.
 - c. START PRIME button pressed a second time prior to priming cycle being completed. The START PRIME button activates and deactivates the timed priming cycle. Pressing the button a second time will immediately stop the priming cycle. Pressing the button a third time will reset the timer and restart the priming cycle from the beginning. Pressing the START PRIME button more than once will result in additional saline being delivered to the device. Always place the device into a holster or over a container to collect the saline that will exit the electrodes as a result of the priming process.
- Keep fingers clear when lowering the pump head to avoid pinching fingers along with the pump segment.

STEP 5: Set RF power

Warnings:

- Use the lowest setting possible to achieve the desired tissue effect to avoid overtreatment, which could result in swelling, fluid, seroma or unintended tissue necrosis.

Set the RF power using the POWER SETTING buttons located on the front panel of the Aquamantys Pump Generator.

- Press the ▲ button to increase the RF power.
- Press the ▼ button to decrease the RF power.

STEP 6: Set saline flow rate

Select the desired saline flow rate setting by pressing the appropriate “Saline Flow Rate” button located on the front panel of the Aquamantys Pump Generator.

Refer to the Aquamantys Pump Generator User Guide for more detailed information on the saline flow rates for each given power setting.

Activating the Aquamantys Disposable Bipolar Device

Warnings:

- Do not activate the Aquamantys disposable bipolar device when the electrodes are not in contact with the tissue to be treated. Activating the device when not in contact with tissue may result in inadvertent tissue damage or user injury due to contact with hot saline.
- Do not wrap electrical cable of device around metal objects. This may induce currents that could lead to shocks, fires, or injuries to the patient or surgical team.

Precautions:

- Confirm that the Aquamantys disposable bipolar device is properly connected to the Aquamantys Pump Generator. Only one device can be connected at any one time.

1. Place the electrodes of the device on the tissue to be treated.
2. Press the activation button on the handpiece of the device to simultaneously activate RF power and saline flow from the device.
3. Release the activation button on the handpiece of the device to shut off both RF power and saline flow from the device.
4. Repeat as necessary by positioning the device over the next area to be treated.
5. Ensure that the electrodes are only in contact with tissue to be coagulated.
6. For optimum performance, the electrodes must be kept free of debris. The electrodes can be cleaned using a sterile moist or dry sponge or sterile brush.
7. If the source of saline needs replenishing, stop using the device and re-prime the saline delivery system and the device using the prime button on the front panel of the generator.
8. If saline flow stops during the electrosurgical procedure, stop using the device and attempt to resume saline flow. Ensure that the pump tubing has been loaded properly into the pump head located on the generator and that the saline bag is not empty. **If unable to resume saline flow, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**

Pressing the activation button on the Aquamantys disposable bipolar device will activate the Aquamantys Pump Generator. The RF Power Activation Indicator will illuminate blue and a continuous RF activation tone will sound to indicate the presence of RF power output. **If the device does not activate when the activation button is depressed, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**

Warnings:

- Do not continuously activate the Aquamantys System for extended periods of time. Extended activation could potentially overheat the pump generator and increase the risk of device malfunction or fire hazard.
- Do not use the Aquamantys disposable bipolar device in an immersed setting (e.g. arthroscopic surgery). Use in an immersed setting could potentially overheat the pump generator and increase the risk of device malfunction or fire hazard.
- Only activate the Aquamantys disposable bipolar device on tissue intended to be treated. Activation over another location may result in hot saline run-off onto unintended tissue, patient, patient drapes, hospital staff, and OR surfaces.
- Use caution to prevent inadvertent activation of the Aquamantys disposable bipolar device during the procedure. Inadvertent activation may result in injury to the patient or surgical team.
- **If the device remains activated when the activation button is released, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**

Changing the RF power setting

Change the RF power setting by using the POWER SETTING buttons located on the front panel of the Aquamantys Pump Generator.

Changing the saline flow rate setting

Precautions:

- Using the low flow rate setting at the high power setting may result in more steam production at the electrodes than with the medium or high flow rate settings, and may result in electrode charring or damage, with reduced hemostatic effectiveness.
- Reduced saline flow can result if one or more of the saline slots or holes in either of the electrodes on the Aquamantys device is clogged by tissue or coagulated blood. If this occurs, clean electrodes with gauze ensuring precautions are taken to avoid inadvertent device activation when cleaning electrodes. Do not attempt to clean electrodes with a scalpel or other surgical instruments. If the above does not correct the problem, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.

The saline flow rate can be adjusted by pressing the appropriate “Saline Flow Rate” button located on the front panel of the Aquamantys Pump Generator.

The Flow Rate Setting Indicator next to the selected flow rate will be illuminated to indicate the current flow rate setting.

Disposing of the Aquamantys Disposable Bipolar Device

Dispose of the Aquamantys disposable bipolar device and the used saline bag according to the procedures for your institution.

Precaution:

- The Aquamantys disposable bipolar device and the saline bag will contain unused saline following use of the device. Take precautions to prevent the unused saline from flowing onto OR surfaces by placing handpiece into waste receptacle prior to opening pump head and removing the pump segment portion of device.

Technical Description

General Description:

Bipolar electrosurgical instrument.
Employs RF energy and saline irrigation for hemostatic sealing and coagulation.
RF energy source: Aquamantys Pump Generator
Generator settings: RF power: 20 – 200 watts
Saline flow rate: Low, Medium, High

General Information:

Sterile, EO
Single-Use, Do Not Reuse
Non-pyrogenic
Does not contain natural rubber latex
Caution: Read Instructions For Use (IFU) before using this device.

Physical Description:

Height: 1.15” (2.92 cm)
Length (without cables): 11” (27.9 cm)
Length of Electrical Cable: approximately 10 feet (3 meters)
Length of Saline Cables: approximately 16 feet (5 meters)
Weight (with cables): 5.0 oz. (143 grams)

Operating Conditions:

Temperature: 50 °F to 86 °F (10 °C to 30 °C)
Humidity: 15% – 85%, non-condensing

Storage Conditions:

Temperature: 50 °F to 104 °F (10 °C to 40 °C)
Humidity: 15% – 85%, non-condensing

Handling Conditions:

Temperature: -20 °F to 122 °F (-29 °C to 50 °C)
Humidity: 15% – 85%, non-condensing

Disposal:

Dispose of used devices as a biohazard.

Limited Express Warranty

IF, PRIOR TO THE PRODUCT EXPIRATION DATE, THE STERILIZED DEVICE IS FOUND TO BE INOPERABLE DURING NORMAL AND PROPER USE IN ACCORDANCE WITH APPLICABLE INSTRUCTIONS, MEDTRONIC ADVANCED ENERGY LLC WILL REPLACE THE PRODUCT AT NO CHARGE. MEDTRONIC ADVANCED ENERGY LLC MAKES NO OTHER WARRANTIES WITH RESPECT TO THE PRODUCT AND EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER. IN NO EVENT SHALL MEDTRONIC ADVANCED ENERGY LLC BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES.

This Limited Express Warranty does not extend to a product where the user has compromised the sterile integrity of the package or the product is used after the product expiration date. A sterile package that has been compromised should not be used.

Customer Service

Please call Medtronic Customer Service at 866.777.9400 if you have any device returns or questions about Medtronic's device.

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